



UNDERSTAND HOW
YOUR BODY WORKS



CE0197



A smart pulse oximeter to spot check your SpO2 and pulse rate



iHealth Air Wireless Pulse Oximeter

iHealth
Connected to my Health

Distributed by

aéras medical

T: 6752 0290

iHealth Air Wireless Pulse Oximeter

The iHealth Pulse Oximeter gives you fast, accurate vital signs readings right from your smartphone or tablet, and is a great and reliable way to spot check your blood oxygen saturation (SpO2) and pulse rate whether you are hiking high at elevations or training for a marathon.

1. Wireless, lightweight, pulse oximeter and easy-to-use and interactive iHealth MyVitals app.
2. Clinically tested and FDA approved.
3. Save data on our secure iHealth Cloud¹ to track results and trends.
4. iHealth MyVitals app supports other iHealth devices for a one-stop overview of all your health-related data.²
5. Compatible with iOS and Android.³



Model PO3M

Color White

Features

PI*/SpO2 ✓

Trend View ✓

Heart Rate ✓

Display LED

Lanyard ✓

Auto Power-OFF ✓

Connectivity

Bluetooth ✓

Cloud ✓

Memory 200

App Supported ✓

iPod/iPad/iPhone ✓

Android Devices ✓

Technical Aspects

Power Battery Rechargeable

Battery Alert ✓

Warranty 1 Year

Dimensions:

27.18mm x 34.54mm x 63.50mm

Weight:

0.04kg

* Perfusion Index.

MyVitals

iHealth



¹ HIPAA compliant.

² iHealth MyVitals app tracks data for iHealth's blood pressure monitors, scales, pulse oximeter and activity trackers.

³ Visit our webpage (aerasmedical.com) for a list of compatible devices.

aéras medical pte ltd

3 Ang Mo Kio Street 62, LINK@AMK #06-20 Singapore 569139

Tel: +65 6752 0290 Fax: +65 6752 6109

Company Registration no.: 201007468Z www.aerasmedical.com