

# MANAGE YOUR BP EFFECTIVELY AND EASILY



**A smart and simple wireless blood pressure monitor that is  
convenient and easy to use by everyone**



**iHealth Feel Wireless Blood Pressure Monitor**

**iHealth**  
Connected to my Health

Distributed by

**aéras medical**

T: 6752 0290

## iHealth Feel Wireless Blood Pressure Monitor

Take control of your health with easy, accurate blood pressure readings anytime, anywhere right from your smartphone or tablet. Our award-winning design and cutting-edge technology makes measuring, storing and tracking your blood pressure virtually effortless.

1. Clinically tested and FDA approved for accurate systolic/diastolic blood pressure and heart rate results every time.
2. Use the iHealth MyVitals app to view results and track trends with intuitive charts/graphs. Measurements evaluated and interpreted using JNC guidelines.<sup>1</sup>
3. One touch data sharing with your doctors and caregivers, with secure iHealth Cloud.<sup>2</sup>
4. iHealth MyVitals app supports other iHealth devices for a one-stop overview of all your health-related data.<sup>3</sup>
5. Wireless, lightweight and award winning design. Comes complete with a travel case.
6. Compatible with iOS and Android.<sup>4</sup>



<b>Model</b>	BP5
<b>Color</b>	White
<b>Features</b>	
<b>Number of User(s)</b>	Unlimited
<b>Pulse Rate</b>	✓
<b>Sensors</b>	Oscillometric
<b>FDA Approved</b>	✓
<b>Connectivity</b>	
<b>Bluetooth</b>	✓
<b>Cloud</b>	✓
<b>Memory</b>	Unlimited
<b>App Supported</b>	✓
<b>iPod/iPad/iPhone<sup>4</sup></b>	✓
<b>Andriod Devices<sup>4</sup></b>	✓
<b>Technical Aspects</b>	
<b>Power Battery</b>	Rechargeable
<b>Low Battery Alert</b>	✓
<b>Warranty</b>	1 Year

### Dimensions:

145mm x 58mm x 30mm

### Weight:

145g

### Cuff circumference:

22-42cm

### Range:

Sys. 60-260 mmHg/ Dia: 40-199 mmHg

Pulse: 40 - 180 bts/min

MyVitals

iHealth



<sup>1</sup> Joint National Committee on Prevention, Detection, Evaluation, and Treatment of High Blood Pressure.

<sup>2</sup> HIPAA compliant.

<sup>3</sup> iHealth MyVitals app tracks data for iHealth's blood pressure monitors, scales, pulse oximeter and activity trackers.

<sup>4</sup> Visit our webpage (aerasmedical.com) for a list of compatible devices.

**aéras medical pte ltd**

3 Ang Mo Kio Street 62, LINK@AMK #06-20 Singapore 569139

Tel: +65 6752 0290 Fax: +65 6752 6109

Company Registration no.: 201007468Z [www.aerasmedical.com](http://www.aerasmedical.com)