


FIRST OF ITS KIND GLUCOMETER



CE0197



A smart blood glucose system that manages diabetes easily



iHealth Wireless Smart Gluco Monitoring System

iHealth
Connected to my Health

Distributed by

aéras medical

T: 6752 0290

iHealth Gluco Wireless Smart Gluco Monitoring System

Doctor trusted, the iHealth Wireless Smart Gluco monitoring system gives you accurate blood sugar readings anytime, anywhere right from your smartphone or tablet. Throw away your log book! This lightweight, portable glucometer works with the iHealth Gluco-smart app to measure and record your readings automatically. Document medication, exercise and diet, and view trends over time for a more complete picture of your health.

1. Small, smart and wireless.
2. The iHealth Gluco-Smart app automatically records and tracks readings. Sync via app to the iHealth Cloud.¹
3. Uses iHealth Blood Glucose Test Strips that work with the iHealth app to track quantity and expiration.
4. Clinically tested/ FDA approved.
5. Meets ISO15197:2013 accuracy standard.
6. Compatible with iOS and Android.²



Model BG5

Color White

Features

Blood Drawn (µL) 0.7

Trend View ✓

Expired Tracking¹ ✓

Screen Display ✓

Test Time 5s

Connectivity

Bluetooth ✓

Cloud ✓

Memory Unlimited

App Supported ✓

Technical Aspects

Display LED

Screen Display ✓

Power Battery Rechargeable

Warranty 1 Year

Dimensions:

14.22mm x 36.57mm x 78.74mm

Weight:

0.03kg

Range Accuracy:

1.1mmol/L - 33.3mmol/L

Operation & Storage:

Operating Temperature: 10 °C - 40 °C

Storage Conditions for Meter and Test Strips:
4° C - 30° C, Humidity 10% - 80% RH

Gluco-Smart



¹ HIPAA compliant.

² Visit our webpage (aerasmedical.com) for a list of compatible devices.

aéras medical pte ltd

3 Ang Mo Kio Street 62, LINK@AMK #06-20 Singapore 569139

Tel: +65 6752 0290 Fax: +65 6752 6109

Company Registration no.: 201007468Z www.aerasmedical.com