

# SMART BP REMOTE MONITORING ANYTIME ANYWHERE



A smart blood pressure monitor easy to use by everyone



**iHealth Track** Connected Arm Blood Pressure Monitor



Distributed by

**aéras medical**

T: 6752 0290

# iHealth®

# Connected to my Health

## iHealth Track Connected Arm Blood Pressure Monitor

Take control of your health with easy, accurate blood pressure readings and heart rate readings with iHealth Track Blood Pressure Monitor. It has been designed for easy use, with or without smartphone. This Blood Pressure Monitor allows you to read and interpret your measurements through its screen that changes colour according to your result.

1. Clinically validated blood pressure by ESH<sup>1</sup> and FDA approved for accurate systolic/diastolic blood pressure and heart rate.
2. Use the iHealth MyVitals app to view results and track trends with intuitive charts/graphs. Measurements evaluated and interpreted using WHO<sup>2</sup> guidelines.
3. One Click: Share all your measurements in your MyVitals app with your doctors and caregivers, with secure iHealth Cloud<sup>3</sup> data storage.
4. iHealth MyVitals app supports other iHealth devices for a one-stop overview of all your health-related data.<sup>4</sup>
5. Easy to use and comfortable.
6. Screen colour changes according to your results.
7. Compatible with iOS and Android.<sup>5</sup>



**Model** KD-550 BT

**Color** White

### Features

**Number of User(s)** Unlimited

**Pulse Rate** ✓

**Sensors** Oscillometric

**Low Battery Alert** ✓

### Connectivity

**Bluetooth** ✓

**Cloud** ✓

**Memory** 99

**App Supported** ✓

**iPod/iPad/iPhone<sup>5</sup>** ✓

**Android Devices<sup>5</sup>** ✓

### Technical Aspects

**ESH<sup>1</sup> Validated** ✓

**Power Battery** 4 x AAA Batteries

**Warranty** 1 Year

### Dimensions:

98mm x 98mm x 46mm

### Weight:

348g (Cuff + Batteries included)

### Cuff circumference:

22-42cm

### Range:

Sys. 60-260 mmHg/ Dia: 40-199 mmHg

Pulse: 40 - 180 bts/min

MyVitals

iHealth

Download on the  
App Store

GET IT ON  
Google Play

Bluetooth™  
SMART

<sup>1</sup>Clinically validated by European Society of Hypertension.

<sup>2</sup>Evaluation of measurements as recommended by the World Health Organisation.

<sup>3</sup>HIPAA compliant.

<sup>4</sup>iHealth MyVitals app tracks data for iHealth's blood pressure monitors, body analysis scales, pulse oximeters and activity trackers.

<sup>5</sup>Visit our webpage (aerasmedical.com) for a list of compatible devices.

**aéras medical pte ltd**

3 Ang Mo Kio Street 62, LINK@AMK #06-20 Singapore 569139

Tel: +65 6752 0290 Fax: +65 6752 6109

Company Registration no.: 201007468Z [www.aerasmedical.com](http://www.aerasmedical.com)