

GLUCOSE METER LIST

CareLink™ iPro™

Medtronic

Some glucose meters supported by this software may not be available in all countries. You will also need a meter cable as supplied by the meter manufacturer. Meter cables are not supplied as part of the iPro2 system. If a patient uses a meter that is not listed here, you can manually enter the BG meter readings into CareLink iPro software.

NOTE: The uploader application in CareLink iPro supports communication with several blood glucose meters from multiple manufacturers, which are listed by Medtronic herein. These devices have been tested and confirmed to work with CareLink iPro. However, please keep in mind: (i) Because Medtronic does not test and confirm every meter in the marketplace, it is possible that other meters (not listed) may also be compatible with CareLink iPro; and (ii) Because Medtronic does not control manufacturers' changes to their meters' hardware, communication protocols, or associated software drivers, Medtronic cannot represent or guarantee continued compatibility of the listed (or non-listed) meters. Therefore, you should regularly confirm the compatibility and functionality of the meter you choose to use with CareLink iPro.

Manufacturer	Blood Glucose Meter Model
Ascensia Diabetes Care	<ul style="list-style-type: none">• CONTOUR®• CONTOUR® NEXT• CONTOUR® NEXT EZ• CONTOUR® NEXT LINK• CONTOUR® NEXT USB• CONTOUR® USB• BREEZE®• BREEZE®2
LifeScan	<ul style="list-style-type: none">• OneTouch® Ultra®• OneTouch® Ultra® 2• OneTouch® UltraLink™• OneTouch® UltraSmart®• OneTouch® UltraMini™
Abbott	<ul style="list-style-type: none">• Precision Xtra®
Roche	<ul style="list-style-type: none">• ACCU-CHEK® Aviva• ACCU-CHEK® Compact Plus• ACCU-CHEK® Nano